

Product datasheet for **TA345402**

ZNF669 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ZNF669 antibody: synthetic peptide directed towards the N terminal of human ZNF669. Synthetic peptide located within the following region: MVSGLRLASRSGEELWLPKPAVARLGPPRHRLRNLRTESPWRSRGSVLFCS
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	52 kDa
Gene Name:	zinc finger protein 669
Database Link:	NP_079080 Entrez Gene 79862 Human Q96BR6
Background:	ZNF669 may be involved in transcriptional regulation.
Synonyms:	FLJ12606
Note:	Immunogen Sequence Homology: Human: 100%
Protein Families:	Transcription Factors



[View online »](#)

Product images:



WB Suggested Anti-ZNF669 Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1: 1562500; Positive
Control: COLO205 cell lysate