


Product datasheet for **TA345408**

ZNF613 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ZNF613 antibody: synthetic peptide directed towards the N terminal of human ZNF613. Synthetic peptide located within the following region: DVMLENYSNLVSVGYQASKPDALFKLEQGEPWTVENEIHSQICPEIKKVN
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	70 kDa
Gene Name:	zinc finger protein 613
Database Link:	NP_001026891 Entrez Gene 79898 Human Q6PF04
Background:	The function of the Anti-ZNF613 gene has not yet been determined.
Synonyms:	FLJ13590
Note:	Immunogen Sequence Homology: Human: 100%; Dog: 85%; Pig: 77%; Horse: 77%
Protein Families:	Transcription Factors


[View online »](#)

Product images:


WB Suggested Anti-ZNF613 Antibody Titration:
 0.2-1 ug/ml; ELISA Titer: 1: 12500; Positive
 Control: Hela cell lysate