

Product datasheet for **TA345502**

zinc finger protein 655 (ZNF655) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ZNF655 antibody: synthetic peptide directed towards the middle region of human ZNF655. Synthetic peptide located within the following region: EGSFSHSSDLILQQEVLTRQKAFDCDVWEKNSSQRAHLVQHQSHTKENS
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	57 kDa
Gene Name:	zinc finger protein 655
Database Link:	NP_612503 Entrez Gene 79027 Human Q8N720
Background:	ZNF655 is a zinc finger protein. The zinc finger proteins are involved in DNA binding and protein-protein interactions. This gene encodes a zinc finger protein. The zinc finger proteins are involved in DNA binding and protein-protein interactions. Multiple alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.
Synonyms:	VIK; VIK-1



[View online »](#)

Note: Immunogen Sequence Homology: Dog: 100%; Human: 100%; Pig: 93%; Rat: 93%; Horse: 93%; Guinea pig: 93%; Mouse: 90%; Bovine: 90%; Rabbit: 90%

Protein Families: Transcription Factors

Product images:



WB Suggested Anti-ZNF655 Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1: 62500; Positive
Control: Human brain