

Product datasheet for **TA345538**


ZNF491 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ZNF491 antibody: synthetic peptide directed towards the N terminal of human ZNF491. Synthetic peptide located within the following region: GERLFESAEGSQCGETFTQVPEDMLNKKTLPGVKSCESGTCGEIFMGYSS
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	51 kDa
Gene Name:	zinc finger protein 491
Database Link:	NP_689569 Entrez Gene 126069 Human Q8N8L2
Background:	The function of ZNF491 has not yet been determined.
Synonyms:	FLJ34791; MGC126639; MGC126641
Note:	Immunogen Sequence Homology: Human: 100%
Protein Families:	Transcription Factors



[View online »](#)

Product images:


WB Suggested Anti-ZNF491 Antibody Titration:
 0.2-1 ug/ml; ELISA Titer: 1: 62500; Positive
 Control: MCF7 cell lysate