

Product datasheet for **TA345547**

ZNF785 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ZNF785 antibody: synthetic peptide directed towards the N terminal of human ZNF785. Synthetic peptide located within the following region: MGPPLAPRPAHVPGEGAPRRTRESRPGAVSFADVAVYFSPEEWECLRPAQ
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	46 kDa
Gene Name:	zinc finger protein 785
Database Link:	NP_689671 Entrez Gene 146540 Human A8K8V0
Background:	ZNF785 may be involved in transcriptional regulation.
Synonyms:	FLJ32130
Note:	Immunogen Sequence Homology: Human: 100%
Protein Families:	Transcription Factors



[View online »](#)

Product images:



WB Suggested Anti-ZNF785 Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1: 62500; Positive
Control: Human Pancreas