

Product datasheet for **TA345580**


ZNF791 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ZNF791 antibody: synthetic peptide directed towards the middle region of human ZNF791. Synthetic peptide located within the following region: CKECGKAFLHSSFQRHTRIHNYEKPLECTQCGKAFLSVSTSLKKHMRMHN
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	67 kDa
Gene Name:	zinc finger protein 791
Database Link:	NP_699189 Entrez Gene 163049 Human Q3KP31
Background:	ZNF791 may be involved in transcriptional regulation.
Synonyms:	FLJ90396; MGC126179; MGC126180
Note:	Immunogen Sequence Homology: Human: 100%; Zebrafish: 82%
Protein Families:	Transcription Factors



[View online »](#)

Product images:

WB Suggested Anti-ZNF791 Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1: 312500; Positive
Control: Human Small Intestine