

Product datasheet for **TA345592**


ZNF57 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ZNF57 antibody: synthetic peptide directed towards the C terminal of human ZNF57. Synthetic peptide located within the following region: LHNHVRMHTGEKPHKCKQCGMSFKWHSSFRNHLRMHTGQKSHECQSYSKA
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification:	Protein A purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	64 kDa
Gene Name:	zinc finger protein 57
Database Link:	NP_775751 Entrez Gene 126295 Human Q68EA5
Background:	The function of ZNF57 remains unknown.
Synonyms:	ZNF424
Note:	Immunogen Sequence Homology: Human: 100%; Horse: 85%; Mouse: 85%; Pig: 83%; Zebrafish: 83%; Rat: 82%
Protein Families:	Transcription Factors


[View online »](#)

Product images:



WB Suggested Anti-ZNF57 Antibody Titration: 5.0 ug/ml; ELISA Titer: 1: 62500; Positive Control: Jurkat cell lysate