

Product datasheet for **TA345727**


GTPBP9 (OLA1) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-GTPBP9 antibody: synthetic peptide directed towards the C terminal of human GTPBP9. Synthetic peptide located within the following region: QAAGKIHTDFEKGFI MAEVMKYEDFKEEGSENAVKAAGKYRQQGRNYIVE
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	44 kDa
Gene Name:	Obg-like ATPase 1
Database Link:	NP_037473 Entrez Gene 29789 Human Q9NTK5
Background:	GTPBP9 belongs to the GTP1/OBG family and the function remains unknown.
Synonyms:	DOC45; GBP45; GTBP9; GTPBP9; PTD004
Note:	Immunogen Sequence Homology: Dog: 100%; Pig: 100%; Rat: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Bovine: 100%; Rabbit: 100%; Zebrafish: 100%; Guinea pig: 100%; Goat: 79%



[View online »](#)

Product images:

WB Suggested Anti-GTPBP9 Antibody Titration:
0.2-1 ug/ml; Positive Control: Jurkat cell lysate