

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# **Product datasheet for TA345855**

# elF4GII (ElF4G3) Rabbit Polyclonal Antibody

# **Product data:**

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
lsotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-EIF4G3 antibody: synthetic peptide directed towards the middle region of human EIF4G3. Synthetic peptide located within the following region: MRGGSSKDLLDNQSQEEQRREMLETVKQLTGGVDVERNSTEAERNKTRES
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. Note that this product is shipped as lyophilized powder to China customers.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	174 kDa
Gene Name:	eukaryotic translation initiation factor 4 gamma 3
Database Link:	<u>NP 003751</u> <u>Entrez Gene 8672 Human</u> <u>O43432</u>
Background:	EIF4G3 is probable component of the protein complex eIF4F, which is involved in the recognition of the mRNA cap, ATP-dependent unwinding of 5'-terminal secondary structure and recruitment of mRNA to the ribosome. It is thought to be a functional homolog of EIF4G1.
Synonyms:	eIF-4G 3; eIF4G 3; eIF4GII




This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

	elF4GII (ElF4G3) Rabbit Polyclonal Antibody – TA345855
Note:	Immunogen Sequence Homology: Dog: 100%; Pig: 100%; Rat: 100%; Horse: 100%; Human:
	100%; Mouse: 100%; Bovine: 100%; Rabbit: 100%; Guinea pig: 100%

Protein Families:Transcription FactorsProtein Pathways:Viral myocarditis

# **Product images:**



WB Suggested Anti-EIF4G3 Antibody Titration: 0.2-1 ug/ml; ELISA Titer: 1: 1562500; Positive Control: Hela cell lysate; EIF4G3 is supported by BioGPS gene expression data to be expressed in HeLa

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US