

## Product datasheet for **TA345874**

### **HNRNPC Rabbit Polyclonal Antibody**

#### **Product data:**

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	WB, IHC
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-HNRNPC antibody: synthetic peptide directed towards the middle region of human HNRNPC. Synthetic peptide located within the following region: ESEGGADDSAEEDLLDDDDNEDRGDDQLELIKDDEKEAEEGEDDRDSAN
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	34 kDa
Gene Name:	heterogeneous nuclear ribonucleoprotein C (C1/C2)
Database Link:	<a href="#">NP_112604</a> <a href="#">Entrez Gene 3183 Human</a> <a href="#">P07910</a>

[View online »](#)

**Background:**

This protein is associated with pre-mRNAs in the nucleus and appears to influence pre-mRNA processing and other aspects of mRNA metabolism and transport. While all of the hnRNPs are present in the nucleus, some seem to shuttle between the nucleus and the cytoplasm. The hnRNP proteins have distinct nucleic acid binding properties. The protein acts as a tetramer and is involved in the assembly of 40S hnRNP particles. This gene belongs to the subfamily of ubiquitously expressed heterogeneous nuclear ribonucleoproteins (hnRNPs). The hnRNPs are RNA binding proteins and they complex with heterogeneous nuclear RNA (hnRNA). These proteins are associated with pre-mRNAs in the nucleus and appear to influence pre-mRNA processing and other aspects of mRNA metabolism and transport. While all of the hnRNPs are present in the nucleus, some seem to shuttle between the nucleus and the cytoplasm. The hnRNP proteins have distinct nucleic acid binding properties. The protein encoded by this gene can act as a tetramer and is involved in the assembly of 40S hnRNP particles. Multiple transcript variants encoding at least two different isoforms have been described for this gene.

**Synonyms:**


C1; C2; HNRNP; HNRPC; SNRPC

**Note:**

Immunogen Sequence Homology: Dog: 100%; Pig: 100%; Rat: 100%; Goat: 100%; Horse: 100%; Human: 100%; Yeast: 100%; Bovine: 100%; Rabbit: 100%; Guinea pig: 100%; Mouse: 93%; Zebrafish: 77%

**Protein Pathways:**

Spliceosome

**Product images:**


IHC Information: Paraffin embedded breast tissue, tested with an antibody Dilution of 5 ug/ml.