

Product datasheet for **TA346131**


cvHSP (HSPB7) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-Hspb7 antibody: synthetic peptide corresponding to a region of Mouse. Synthetic peptide located within the following region: TFAHKCQLPEDVDPTSVTLSALREDGSLTIRARRNPHEHVQQTFRTEIKI
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	19 kDa
Gene Name:	heat shock protein family B (small) member 7
Database Link:	NP_055239 Entrez Gene 29818 Mouse Q9UBY9
Background:	The function of Hspb7 remains unknown.
Synonyms:	cvHSP
Note:	Immunogen Sequence Homology: Dog: 100%; Pig: 100%; Rat: 100%; Human: 100%; Mouse: 100%; Bovine: 100%; Rabbit: 100%; Guinea pig: 100%; Horse: 93%; Zebrafish: 83%



[View online »](#)

Product images:

WB Suggested Anti-Hspb7 Antibody Titration: 0.2-1 ug/ml; ELISA Titer: 1: 62500; Positive Control: Mouse Heart