

Product datasheet for TA346281

CD40 Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: WB

Recommended Dilution: WB

Reactivity: Human

Host: Rabbit

Isotype: IgG

Clonality: Polyclonal

Immunogen: The immunogen for anti-CD40 antibody: synthetic peptide directed towards the N terminal of

human CD40. Synthetic peptide located within the following region: SQCCSLCQPGQKLVSDCTEFTETECLPCGESEFLDTWNRETHCHQHKYCD

Formulation: Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2%

sucrose.

Note that this product is shipped as lyophilized powder to China customers.

Purification: Affinity Purified

Conjugation: Unconjugated

Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 20 kDa

Gene Name: CD40 molecule

Database Link: NP 690593

Entrez Gene 958 Human

P25942



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



Background: CD40 is the receptor for TNFSF5/CD40LG. Defects in CD40 are the cause of hyper-IgM

immunodeficiency type 3 (HIGM3). The protein encoded by this gene is a member of the TNF-receptor superfamily. This receptor has been found to be essential in mediating a broad variety of immune and inflammatory responses including T cell-dependent immunoglobulin class switching, memory B cell development, and germinal center formation. AT-hook transcription factor AKNA is reported to coordinately regulate the expression of this receptor and its ligand, which may be important for homotypic cell interactions. Adaptor protein TNFR2 interacts with this receptor and serves as a mediator of the signal transduction. The interaction of this receptor and its ligand is found to be necessary for amyloid-beta-induced microglial activation, and thus is thought to be an early event in Alzheimer disease pathogenesis. Two alternatively spliced transcript variants of this gene encoding distinct isoforms have been reported.

Synonyms: Bp50; CDW40; p50; TNFRSF5

Note: Immunogen Sequence Homology: Human: 100%; Rabbit: 100%; Pig: 91%; Guinea pig: 91%;

Horse: 83%


Protein Families: Druggable Genome, Secreted Protein, Transmembrane

Protein Pathways: Allograft rejection, Asthma, Autoimmune thyroid disease, Cell adhesion molecules (CAMs),

Cytokine-cytokine receptor interaction, Primary immunodeficiency, Systemic lupus

erythematosus, Toll-like receptor signaling pathway, Viral myocarditis

Product images:



WB Suggested Anti-CD40 Antibody Titration: 0.2-1 ug/ml; ELISA Titer: 1: 62500; Positive Control:

Human Thymus