

Product datasheet for TA346469

PPP2R3B Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: WB

Recommended Dilution: WB

Reactivity: Human

Host: Rabbit

Isotype: IgG

Clonality: Polyclonal

Immunogen: The immunogen for anti-PPP2R3B antibody: synthetic peptide directed towards the C

terminal of human PPP2R3B. Synthetic peptide located within the following region:

ELSDWEKYAAEEYDILVAEETAGEPWEDGFEAELSPVEQKLSALRSPLAQ

Formulation: Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2%

sucrose.

Note that this product is shipped as lyophilized powder to China customers.

Purification: Affinity Purified

Conjugation: Unconjugated

Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 65 kDa

Gene Name: protein phosphatase 2 regulatory subunit B"beta

Database Link: NP 037371

Entrez Gene 28227 Human

Q9Y5P8



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



Background:


Protein phosphatase 2 (formerly named type 2A) is one of the four major Ser/Thr phosphatases and is implicated in the negative control of cell growth and division. Protein phosphatase 2 holoenzymes are heterotrimeric proteins composed of a structural subunit A, a catalytic subunit C, and a regulatory subunit B. The regulatory subunit is encoded by a diverse set of genes that have been grouped into the B/PR55, B'/PR61, and B"/PR72 families. These different regulatory subunits confer distinct enzymatic specificities and intracellular localizations to the holozenzyme. PPP2R3B belongs to the B" family. The B" family has been further divided into subfamilies. PPP2R3B belongs to the beta subfamily of regulatory subunit B".Protein phosphatase 2 (formerly named type 2A) is one of the four major Ser/Thr phosphatases and is implicated in the negative control of cell growth and division. Protein phosphatase 2 holoenzymes are heterotrimeric proteins composed of a structural subunit A, a catalytic subunit C, and a regulatory subunit B. The regulatory subunit is encoded by a diverse set of genes that have been grouped into the B/PR55, B'/PR61, and B"/PR72 families. These different regulatory subunits confer distinct enzymatic specificities and intracellular localizations to the holozenzyme. The product of this gene belongs to the B" family. The B" family has been further divided into subfamilies. The product of this gene belongs to the beta subfamily of regulatory subunit B". Alternative splicing results in multiple transcript variants encoding different isoforms.

Synonyms: NYREN8; PPP2R3L; PPP2R3LY; PR48

Note: Immunogen Sequence Homology: Human: 100%; Dog: 93%; Bovine: 86%; Rat: 75%

Protein Families: Druggable Genome, Phosphatase

Product images:



WB Suggested Anti-PPP2R3B Antibody Titration: 0.2-1 ug/ml; Positive Control: HepG2 cell lysate