

Product datasheet for **TA346627**

TEX10 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-TEX10 antibody is: synthetic peptide directed towards the middle region of Human TEX10. Synthetic peptide located within the following region: RLTSQQWRLKVLVRLSKFLQALADGSSRLRESEGLQEQKENPHATSNSIF
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	47 kDa
Gene Name:	testis expressed 10
Database Link:	NP_001155056 Entrez Gene 54881 Human Q9NXF1
Background:	The function of this protein remains unknown.
Synonyms:	bA208F1.2; Ipi1
Note:	Immunogen Sequence Homology: Dog: 100%; Pig: 100%; Rat: 100%; Horse: 100%; Human: 100%; Bovine: 100%; Rabbit: 100%; Guinea pig: 100%; Mouse: 93%



[View online »](#)

Product images:



Host: Rabbit; Target Name: TEX10; Sample
Tissue: Fetal Brain lysates; Antibody Dilution: 1.0
ug/ml