

## Product datasheet for **TA346787**


### Ftsj3 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-Ftsj3 antibody: synthetic peptide corresponding to a region of Rat. Synthetic peptide located within the following region: QVTYVVAKKGVGRKVRPAGVRGHFKVVDSRMKKDQRAQRKEQKRNHRRK
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	100 kDa
Gene Name:	Ftsj homolog 3
Database Link:	<a href="#">NP_001012014</a> <a href="#">Entrez Gene 303608 Rat</a> <a href="#">Q5RJT2</a>
Background:	The function of Ftsj3 remains unknown.
Synonyms:	EPCS3; FLJ20062
Note:	Immunogen Sequence Homology: Dog: 100%; Pig: 100%; Rat: 100%; Human: 100%; Mouse: 100%; Bovine: 100%; Guinea pig: 100%; Horse: 93%; Rabbit: 93%



[View online »](#)

**Product images:**

WB Suggested Anti-Ftsj3 Antibody Titration: 0.2-1 ug/ml; ELISA Titer: 1: 12500; Positive Control: Rat Heart