

Product datasheet for **TA346853**

C3orf39 (POMGNT2) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-C3orf39 antibody: synthetic peptide directed towards the N terminal of human C3orf39. Synthetic peptide located within the following region: MHLSAVFNALLVSVLA AVLWKHVRLREHAATLEEELALSRQATEPAPALR
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	66 kDa
Gene Name:	protein O-linked mannose N-acetylglucosaminyltransferase 2 (beta 1,4-)
Database Link:	NP_116195 Entrez Gene 84892 Human Q8NAT1
Background:	This gene encodes a protein with glycosyltransferase activity although its function is not currently known. [provided by RefSeq, Sep 2012]
Synonyms:	AGO61; C3orf39; GTDC2; MDDGA8
Note:	Immunogen Sequence Homology: Rat: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Rabbit: 100%; Dog: 93%; Pig: 93%; Bovine: 93%; Guinea pig: 93%



[View online »](#)

Product images:



WB Suggested Anti-C3orf39 Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1: 312500; Positive
Control: SH-SY5Y cell lysate