

Product datasheet for **TA346897S**

BRAP Mouse Monoclonal Antibody [Clone ID: 1E7-C9-D10-B10]

Product data:

Product Type:	Primary Antibodies
Clone Name:	1E7-C9-D10-B10
Applications:	WB
Recommended Dilution:	WB: 1:1000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	The immunogen for BRAP antibody: purified recombinant human BRAP protein fragments expressed in E.coli.
Formulation:	Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	67 kDa
Gene Name:	BRCA1 associated protein
Database Link:	NP_006759 Entrez Gene 8315 Human Q7Z569
Background:	The protein encoded by this gene was identified by its ability to bind to the nuclear localization signal of BRCA1 and other proteins. It is a cytoplasmic protein which may regulate nuclear targeting by retaining proteins with a nuclear localization signal in the cytoplasm. [provided by RefSeq, Jul 2008]
Synonyms:	BRAP2; IMP; RNF52
Protein Families:	Druggable Genome

[View online »](#)