

## Product datasheet for **TA346949**

### FOXP3 Mouse Monoclonal Antibody [Clone ID: 7H9-D6-A10]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	7H9-D6-A10
Applications:	IF, WB
Recommended Dilution:	WB: 1:300, IF: 1:100
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	The immunogen for FOXP3 antibody: purified recombinant human FOXP3 protein fragments expressed in E.coli.
Formulation:	Ascites containing with 0.02% sodium azide and 50% glycerol.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	47 kDa
Gene Name:	forkhead box P3
Database Link:	<a href="#">NP_054728</a> <a href="#">Entrez Gene 50943 Human</a> <a href="#">Q9BZS1</a>
Background:	The protein encoded by this gene is a member of the forkhead/winged-helix family of transcriptional regulators. Defects in this gene are the cause of immunodeficiency polyendocrinopathy, enteropathy, X-linked syndrome (IPEX), also known as X-linked autoimmunity-immunodeficiency syndrome. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2008]
Synonyms:	AIID; DIETER; IPEX; JM2; PIDX; XPID
Protein Families:	Transcription Factors



[View online »](#)