

Product datasheet for **TA346953S**

GATA6 Mouse Monoclonal Antibody [Clone ID: 3H3-H6-H11]

Product data:

Product Type:	Primary Antibodies
Clone Name:	3H3-H6-H11
Applications:	WB
Recommended Dilution:	WB: 1:200
Reactivity:	Human
Host:	Mouse
Clonality:	Monoclonal
Immunogen:	The immunogen for GATA6 antibody: purified recombinant human GATA6 protein fragments expressed in E.coli.
Formulation:	ascites
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	60 kDa
Gene Name:	GATA binding protein 6
Database Link:	NP_005248 Entrez Gene 2627 Human Q92908
Background:	Thought to be important for regulating terminal differentiation and/or proliferation.
Synonyms:	GATA-binding protein 6; GATA binding protein 6; transcription factor GATA-6
Protein Families:	Embryonic stem cells, ES Cell Differentiation/IPS, Transcription Factors



[View online »](#)