

## Product datasheet for **TA346956**

### TFIIB (GTF2B) Mouse Monoclonal Antibody [Clone ID: 2F6-A3]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	2F6-A3
Applications:	WB
Recommended Dilution:	WB: 1:500
Reactivity:	Human, Monkey, Mouse
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Immunogen:	The immunogen for GTF2B antibody: purified recombinant human TFIIB protein fragments expressed in E.coli.
Formulation:	ascites
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	34 kDa
Gene Name:	general transcription factor IIB
Database Link:	<a href="#">NP_001505</a> <a href="#">Entrez Gene 229906 Mouse</a> <a href="#">Entrez Gene 2959 Human</a> <a href="#">Q00403</a>
Background:	GTF2B is one of the ubiquitous factors required for transcription initiation by RNA polymerase II. The protein localizes to the nucleus where it forms a complex (the DAB complex) with transcription factors IID and IIA. Transcription factor IIB serves as a bridge between IID, the factor which initially recognizes the promoter sequence, and RNA polymerase II.
Synonyms:	TF2B; TFIIB
Protein Families:	Druggable Genome, Transcription Factors
Protein Pathways:	Basal transcription factors



[View online »](#)