

Product datasheet for TA347021

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

PPP1A (PPP1CA) Mouse Monoclonal Antibody [Clone ID: 1C11-C10-H5-E9]

Product data:

Product Type: Primary Antibodies

Clone Name: 1C11-C10-H5-E9

Applications: IF, IHC, WB

Recommended Dilution: WB: 1:1000, IF: 1:300

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: The immunogen for PPP1CA antibody: purified recombinant human PPP1A protein fragments

expressed in E.coli.

Formulation: Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with

0.02% sodium azide, 50% glycerol

Purification: Affinity purified Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 38 kDa

Gene Name: protein phosphatase 1 catalytic subunit alpha

Database Link: NP 001008709

Entrez Gene 5499 Human

P62136





PPP1A (PPP1CA) Mouse Monoclonal Antibody [Clone ID: 1C11-C10-H5-E9] - TA347021

Background: The protein encoded by this gene is one of the three catalytic subunits of protein

phosphatase 1 (PP1). PP1 is a serine/threonine specific protein phosphatase known to be involved in the regulation of a variety of cellular processes, such as cell division, glycogen metabolism, muscle contractility, protein synthesis, and HIV-1 viral transcription. Increased PP1 activity has been observed in the end stage of heart failure. Studies in both human and mice suggest that PP1 is an important regulator of cardiac function. Mouse studies also suggest that PP1 functions as a suppressor of learning and memory. Three alternatively spliced transcript variants encoding different isoforms have been found for this gene.

[provided by RefSeq, Jul 2008]

Synonyms: PP-1A; PP1A; PP1alpha; PPP1A

Protein Families: Druggable Genome, Phosphatase

Protein Pathways: Focal adhesion, Insulin signaling pathway, Long-term potentiation, Oocyte meiosis, Regulation

of actin cytoskeleton, Vascular smooth muscle contraction