

## Product datasheet for **TA347025S**

### PPP4C Mouse Monoclonal Antibody [Clone ID: 2F11-D10-G4]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	2F11-D10-G4
Applications:	IF, IHC, WB
Recommended Dilution:	WB: 1:200, IF: 1:200
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	The immunogen for PPP4C antibody: purified recombinant human Protein Phosphatase 4C protein fragments expressed in E.coli.
Formulation:	Purified mouse monoclonal in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	34 kDa
Gene Name:	protein phosphatase 4 catalytic subunit
Database Link:	<a href="#">NP_002711</a> <a href="#">Entrez Gene 5531 Human</a> <a href="#">P60510</a>
Background:	Protein phosphatase 4C may be involved in microtubule organization. It binds 1 iron ion and 1 manganese ion per subunit. PP4 consists of a catalytic subunit PPP4C and a regulatory subunit. PPP4R1 and belongs to the PPP phosphatase family, PP X subfamily.
Synonyms:	PP-X; PP4; PP4C; PPH3; PPP4; PPX
Protein Families:	Druggable Genome, Phosphatase



[View online »](#)