

## Product datasheet for **TA347028S**

### DNA PKcs (PRKDC) Mouse Monoclonal Antibody [Clone ID: 6D1-C11-F10]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	6D1-C11-F10
Applications:	IF, IHC, IP, WB
Recommended Dilution:	WB: 1:1000, IF: 1:100
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	The immunogen for PRKDC antibody: purified recombinant human DNA-PKcs protein fragments expressed in E.coli
Formulation:	Purified mouse monoclonal antibody in PBS(pH 7.4)containing with 0.02% sodium azide,0.1%BSA & 50% glycerol.
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	450 kDa
Gene Name:	protein kinase, DNA-activated, catalytic polypeptide
Database Link:	<a href="#">NP_008835</a> <a href="#">Entrez Gene 5591 Human</a> <a href="#">P78527</a>
Background:	This gene encodes the catalytic subunit of the DNA-dependent protein kinase(DNA-PK).It functions with the Ku70/Ku80 heterodimer protein in DNA double strand break repair and recombination.The protein encoded is a member of the PI3/PI4-kinase family. [provided by RefSeq,Jul 2010]
Synonyms:	DNA-PKcs; DNAPK; DNPk1; HYRC; HYRC1; IMD26; p350; XRCC7
Protein Families:	Druggable Genome, Protein Kinase



[View online »](#)

Protein Pathways: Cell cycle, Non-homologous end-joining