

Product datasheet for **TA347055M**

SMYD4 Mouse Monoclonal Antibody [Clone ID: 7C8-B11-F10]

Product data:

Product Type:	Primary Antibodies
Clone Name:	7C8-B11-F10
Applications:	WB
Recommended Dilution:	WB: 1:200
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	The immunogen for SMYD4 antibody: purified recombinant human SMYD4 protein fragments expressed in E.coli.
Formulation:	ascites
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	89 kDa
Gene Name:	SET and MYND domain containing 4
Database Link:	NP_443160 Entrez Gene 114826 Human Q8IYR2
Synonyms:	ZMYND21



[View online »](#)