

Product datasheet for **TA347063**

TIF1 gamma (TRIM33) Mouse Monoclonal Antibody [Clone ID: 3H7-C8-B9]

Product data:

Product Type:	Primary Antibodies
Clone Name:	3H7-C8-B9
Applications:	WB
Recommended Dilution:	WB: 1:1000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	The immunogen for TRIM33 antibody: purified recombinant human TIF1 gamma / TRIM33 protein fragments expressed in E.coli.
Formulation:	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.02% sodium azide, 50% glycerol
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	140 kDa
Gene Name:	tripartite motif containing 33
Database Link:	NP_056990 Entrez Gene 51592 Human Q9UPN9
Background:	The protein encoded by this gene is thought to be a transcriptional corepressor. However, molecules that interact with this protein have not yet been identified. The protein is a member of the tripartite motif family. This motif includes three zinc-binding domains, a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region. Three alternatively spliced transcript variants for this gene have been described, however, the full-length nature of one variant has not been determined.
Synonyms:	ECTO; PTC7; RFG7; TF1G; TIF1G; TIF1GAMMA; TIFGAMMA



[View online »](#)

Protein Families: Druggable Genome, ES Cell Differentiation/IPS, Protein Kinase, Transcription Factors