

Product datasheet for TA347331

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

14-3-3 zeta (YWHAZ) Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: WB

Reactivity: WB: 1:500-1:3000 Human, Mouse, Rat

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

Immunogen: The immunogen for anti-14-3-3 beta/? Antibody: A synthesized peptide derived from human

14-3-3 beta/?

Formulation: Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50%

glycerol. Store at -20?. Stable for 12 months from date of receipt

Concentration: lot specific

Purification: Affinity-chromatography

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 28 kDa

Gene Name: tyrosine 3-monooxygenase/tryptophan 5-monooxygenase activation protein zeta

Database Link: NP 003397

Entrez Gene 22631 MouseEntrez Gene 25578 RatEntrez Gene 7534 Human

P63104

Background: 14-3-3 zeta a protein of the 14-3-3 family of proteins which mediate signal transduction by

binding to phosphoserine-containing proteins. A multifunctional regulator of the cell signaling

processes. Phosphorylation apparently disrupts homodimerization.

Synonyms: 14-3-3-zeta; HEL-S-3; HEL4; KCIP-1; YWHAD

Note: 14-3-3 beta/? antibody detects endogenous levels of total 14-3-3 beta/?






Protein Pathways:

Cell cycle, Neurotrophin signaling pathway, Oocyte meiosis, Pathogenic Escherichia coli infection

Product images:



Western blot analysis on RAW264.7 cell lysate using 14-3-3 beta/? Antibody