

Product datasheet for **TA347452**

Parathyroid Hormone (PTH) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500-1:3000
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-Parathyroid Hormone Antibody: A synthesized peptide derived from human Parathyroid Hormone
Formulation:	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Store at -20?. Stable for 12 months from date of receipt
Concentration:	lot specific
Purification:	Affinity-chromatography
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	16 kDa
Gene Name:	parathyroid hormone
Database Link:	NP_000306 Entrez Gene 5741 Human P01270
Background:	PTH PTH elevates calcium level by dissolving the salts in bone and preventing their renal excretion. Belongs to the parathyroid hormone family.
Synonyms:	PTH1
Note:	Parathyroid Hormone antibody detects endogenous levels of total Parathyroid Hormone
Protein Families:	Druggable Genome, Secreted Protein



[View online »](#)

Product images:



Western blot analysis on A549 cell lysate using Parathyroid Hormone Antibody