

Product datasheet for **TA347829**

FOXP4 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500-1:1000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.Store at -20 ?.Stable for 12 months from date of receipt
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	74kd
Gene Name:	forkhead box P4
Database Link:	NP_612466 Entrez Gene 74123 MouseEntrez Gene 363185 RatEntrez Gene 116113 Human Q8IVH2
Background:	Transcriptional repressor that represses lung-specific expression.
Synonyms:	hFKHLA
Note:	The antibody detects endogenous level of total FOXP4 protein.
Protein Families:	Transcription Factors



[View online »](#)

Product images:



Western blot analysis of extract from JK and HeLa cells using FOXP4