

Product datasheet for **TA348019**

HNRNPD Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500-1:2000
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-hnRPD Antibody: A synthesized peptide derived from human hnRPD
Formulation:	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.Store at -20 ?.Stable for 12 months from date of receipt
Concentration:	lot specific
Purification:	Immunogen affinity purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	38 kDa
Gene Name:	heterogeneous nuclear ribonucleoprotein D
Database Link:	NP_112737 Entrez Gene 11991 Mouse Entrez Gene 3184 Human Q14103
Background:	Binds with high affinity to RNA molecules that contain AU-rich elements (AREs) found within the 3' UTR of many proto-oncogenes and cytokine mRNAs. Also binds to double- and single-stranded DNA sequences in a specific manner and functions a transcription factor.
Synonyms:	AUF1; AUF1A; hnRNPD0; HNRPD; P37
Note:	hnRPD Antibody detects endogenous levels of total hnRPD
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)

Product images:



Western blot analysis of hnRPD expression in 293 cell extract