

Product datasheet for **TA349652**

FXYD1 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 500-2000 WB positive control: Mouse heart tissue IHC: 50-200 Positive control: Human tonsil Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human FXYD1
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	10 kDa
Gene Name:	FXYD domain containing ion transport regulator 1
Database Link:	NP_068702 Entrez Gene 58971 Rat Entrez Gene 5348 Human O00168



[View online »](#)

Background:


This gene encodes a member of a family of small membrane proteins that share a 35-amino acid signature sequence domain, beginning with the sequence PFXYD and containing 7 invariant and 6 highly conserved amino acids. The approved human gene nomenclature for the family is FXYD-domain containing ion transport regulator. Mouse FXYD5 has been termed RIC (Related to Ion Channel). FXYD2, also known as the gamma subunit of the Na,K-ATPase, regulates the properties of that enzyme. FXYD1 (phospholemman), FXYD2 (gamma), FXYD3 (MAT-8), FXYD4 (CHIF), and FXYD5 (RIC) have been shown to induce channel activity in experimental expression systems.

Synonyms:

PLM

Protein Families:

Ion Channels: Other, Transmembrane

Product images:

Gel: 12%SDS-PAGE


Lysate: 40 µg

Lane: Mouse heart tissue


Primary antibody: TA349652 (FXYD1 Antibody) at dilution 1/1290

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution


Exposure time: 1 second




Immunohistochemistry of paraffin-embedded Human tonsil tissue using TA349652 (FXYD1 Antibody) at dilution 1/60 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human tonsil tissue using TA349652 (FXYD1 Antibody) at dilution 1/60, treated with fusion protein. (Original magnification: x200)



Immunohistochemistry of paraffin-embedded Human colon cancer tissue using TA349652 (FXYD1 Antibody) at dilution 1/60 (Original magnification: x200)



Immunohistochemistry of paraffin-embedded Human colon cancer tissue using TA349652 (FXYD1 Antibody) at dilution 1/60, treated with fusion protein. (Original magnification: x200)