

Product datasheet for TA349654

CERK Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: IHC

Recommended Dilution: IHC: 25-100

Positive control: Human tonsil Predicted cell location: Cytoplasm

Reactivity: Human, Mouse

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

Immunogen: Fusion protein of human CERK

Formulation: pH7.4 PBS, 0.05% NaN3, 40% Glyceroln

Concentration: lot specific

Purification: Antigen affinity purification

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Gene Name: ceramide kinase

Database Link: NP 073603

Entrez Gene 64781 Human

Q8TCT0

Background: CERK converts ceramide to ceramide 1-phosphate (C1P), a sphingolipid metabolite. Both

CERK and C1P have been implicated in various cellular processes, including proliferation, apoptosis, phagocytosis, and inflammation. Catalyzes specifically the phosphorylation of ceramide to form ceramide 1-phosphate. Acts efficiently on natural and analog ceramides (C6, C8, C16 ceramides, and C8-dihydroceramide), to a lesser extent on C2-ceramide and C6-

dihydroceramide, but not on other lipids, such as various sphingosines. Binds

phosphoinositides.

Synonyms: dA59H18.2; dA59H18.3; hCERK; LK4



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn


Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com

CERK Rabbit Polyclonal Antibody - TA349654


Protein Families: Druggable Genome

Protein Pathways: Metabolic pathways, Sphingolipid metabolism

Product images:




Immunohistochemistry of paraffin-embedded Human tonsil tissue using TA349654 (CERK Antibody) at dilution 1/50 (Original magnification: ×200)




Immunohistochemistry of paraffin-embedded Human tonsil tissue using TA349654 (CERK Antibody) at dilution 1/50, treated with fusion protein. (Original magnification: ×200)





Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using TA349654 (CERK Antibody) at dilution 1/50 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using TA349654 (CERK Antibody) at dilution 1/50, treated with fusion protein. (Original magnification: ×200)