

Product datasheet for TA350108

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

gamma Catenin (JUP) Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: IHC, WB

Recommended Dilution: WB: 1000-5000

WB positive control: Mouse heart tissue and Hela cells, HUVEC cells and mouse skin tissue

IHC: 50-200

Positive control: Human cervical cancer Predicted cell location: Cytoplasm

Reactivity: Human, Mouse, Rat

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

Immunogen: Fusion protein of human JUP

Formulation: pH7.4 PBS, 0.05% NaN3, 40% Glyceroln

Concentration: lot specific

Purification: Antigen affinity purification

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 82 kDa

Gene Name: junction plakoglobin

Database Link: NP 002221

Entrez Gene 16480 MouseEntrez Gene 81679 RatEntrez Gene 3728 Human

P14923





Background: This gene encodes a major cytoplasmic protein which is the only known constituent common

to submembranous plaques of both desmosomes and intermediate junctions. This protein forms distinct complexes with cadherins and desmosomal cadherins and is a member of the catenin family since it contains a distinct repeating amino acid motif called the armadillo repeat. Mutation in this gene has been associated with Naxos disease. Alternative splicing

occurs in this gene; however, not all transcripts have been fully described.


Synonyms: ARVD12; CTNNG; DP3; DPIII; PDGB; PKGB

Protein Families: Druggable Genome

Protein Pathways: Acute myeloid leukemia, Arrhythmogenic right ventricular cardiomyopathy (ARVC), Pathways

in cancer

Product images:



Gel: 6%SDS-PAGE Lysate: 40 μg

Lane 1-4: Mouse heart tissue

Hela cells HUVEC cells mouse skin tissue


Primary antibody: TA350108 (JUP Antibody) at

dilution 1/1150

Secondary antibody: Goat anti rabbit IgG at


1/8000 dilution

Exposure time: 10 seconds




Immunohistochemistry of paraffin-embedded Human cervical cancer tissue using TA350108 (JUP Antibody) at dilution 1/50 (Original magnification: ×200)






Immunohistochemistry of paraffin-embedded Human cervical cancer tissue using TA350108 (JUP Antibody) at dilution 1/50, treated with fusion protein. (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human lung cancer tissue using TA350108 (JUP Antibody) at dilution 1/50 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human lung cancer tissue using TA350108 (JUP Antibody) at dilution 1/50, treated with fusion protein. (Original magnification: ×200)