

Product datasheet for TA350206

NDUFA1 Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: IHC, WB

Recommended Dilution: WB: 500-2000

WB positive control: 293T cell, Mouse brain tissue, Mouse heart tissue, Rat skeletal muscle

tissue lysates IHC: 100-300

Positive control: Human breast cancer Predicted cell location: Cytoplasm

Reactivity: Human, Mouse, Rat

Rabbit Host: Isotype: lgG

Clonality: Polyclonal

Immunogen: Fusion protein of human NDUFA1

Formulation: pH7.4 PBS, 0.05% NaN3, 40% Glyceroln

Concentration: lot specific

Purification: Antigen affinity purification

Conjugation: Unconjugated

Store at -20°C as received. Storage:

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 8 kDa

Gene Name: NADH:ubiquinone oxidoreductase subunit A1

Database Link: NP 004532

Entrez Gene 4694 Human

O15239



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



Background:

The human NDUFA1 gene codes for an essential component of complex I of the respiratory chain, which transfers electrons from NADH to ubiquinone. It has been noted that the N-terminal hydrophobic domain has the potential to be folded into an alpha-helix spanning the inner mitochondrial membrane with a C-terminal hydrophilic domain interacting with globular subunits of complex I. The highly conserved two-domain structure suggests that this feature is critical for the protein function and might act as an anchor for the NADH:ubiquinone oxidoreductase complex at the inner mitochondrial membrane.


Synonyms: CI-MWFE; MWFE; ZNF183

Protein Families: Transmembrane

Protein Pathways: Alzheimer's disease, Huntington's disease, Metabolic pathways, Oxidative phosphorylation,

Parkinson's disease

Product images:



Gel: 12%SDS-PAGE Lysate: 40 µg Lane 1-4: 293T cell Mouse brain tissue Mouse heart tissue

Rat skeletal muscle tissue lysates


Primary antibody: TA350206 (NDUFA1 Antibody)

at dilution 1/800

Secondary antibody: Goat anti rabbit IgG at


1/5000 dilution

Exposure time: 40 seconds




Immunohistochemistry of paraffin-embedded Human breast cancer tissue using TA350206 (NDUFA1 Antibody) at dilution 1/100 (Original magnification: ×200)






Immunohistochemistry of paraffin-embedded Human breast cancer tissue using TA350206 (NDUFA1 Antibody) at dilution 1/100, treated with fusion protein. (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human gastric cancer tissue using TA350206 (NDUFA1 Antibody) at dilution 1/100 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human gastric cancer tissue using TA350206 (NDUFA1 Antibody) at dilution 1/100, treated with fusion protein. (Original magnification: ×200)