

## Product datasheet for **TA350250S**

### ODF2 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 500-2000 WB positive control: Human testis tissue and K562 cells IHC: 25-100 Positive control: Human colon cancer Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human ODF2
Formulation:	pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	95 kDa
Gene Name:	outer dense fiber of sperm tails 2
Database Link:	<a href="#">NP_002531</a> <a href="#">Entrez Gene 18286 Mouse</a> <a href="#">Entrez Gene 29479 Rat</a> <a href="#">Entrez Gene 4957 Human</a> <a href="#">Q5BJF6</a>


**Background:** The outer dense fibers are cytoskeletal structures that surround the axoneme in the middle piece and principal piece of the sperm tail. The fibers function in maintaining the elastic structure and recoil of the sperm tail as well as in protecting the tail from shear forces during epididymal transport and ejaculation. Defects in the outer dense fibers lead to abnormal sperm morphology and infertility. This gene encodes one of the major outer dense fiber proteins. Alternative splicing results in multiple transcript variants.




[View online »](#)

Synonyms: 1; 2; CT134; ODF2; ODF84


**Product images:**



Gel: 6%SDS-PAGE  
 Lysate: 40 µg  
 Lane 1-2: Human testis tissue  
 K562 cells  
 Primary antibody: [TA350250] (ODF2 Antibody) at dilution 1/650  
 Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution  
 Exposure time: 30 seconds



Immunohistochemistry of paraffin-embedded Human colon cancer tissue using [TA350250] (ODF2 Antibody) at dilution 1/30 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human colon cancer tissue using [TA350250] (ODF2 Antibody) at dilution 1/30, treated with fusion protein. (Original magnification: ×200)