

Product datasheet for **TA350272S**

PDLIM1 Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 500-2000 WB positive control: A172 cells IHC: 25-100 Positive control: Human colon cancer Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human PDLIM1
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	36 kDa
Gene Name:	PDZ and LIM domain 1
Database Link:	NP_066272 Entrez Gene 54132 Mouse Entrez Gene 54133 Rat Entrez Gene 9124 Human O00151

Background: This gene encodes a member of the enigma protein family. The protein contains two protein interacting domains, a PDZ domain at the amino terminal end and one to three LIM domains at the carboxyl terminal. It is a cytoplasmic protein associated with the cytoskeleton. The protein may function as an adapter to bring other LIM-interacting proteins to the cytoskeleton. Pseudogenes associated with this gene are located on chromosomes 3, 14 and 17.



[View online »](#)

Synonyms: CLIM1; CLP-36; CLP36; hCLIM1; HEL-S-112


Product images:



Gel: 8%SDS-PAGE
Lysate: 40 µg
Lane: A172 cells
Primary antibody: [TA350272] (PDLIM1 Antibody) at dilution 1/400
Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution
Exposure time: 10 seconds



Immunohistochemistry of paraffin-embedded Human colon cancer tissue using [TA350272] (PDLIM1 Antibody) at dilution 1/30 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human colon cancer tissue using [TA350272] (PDLIM1 Antibody) at dilution 1/30, treated with fusion protein. (Original magnification: ×200)