

Product datasheet for **TA350319**

RAB27A Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 500-2000 WB positive control: K562 cells IHC: 25-100 Positive control: Human liver cancer Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human RAB27A
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	25 kDa
Gene Name:	RAB27A, member RAS oncogene family
Database Link:	NP_899057 Entrez Gene 11891 Mouse Entrez Gene 50645 Rat Entrez Gene 5873 Human P51159

Background: The protein encoded by this gene belongs to the small GTPase superfamily, Rab family. The protein is membrane-bound and may be involved in protein transport and small GTPase mediated signal transduction. Mutations in this gene are associated with Griscelli syndrome type 2. Alternative splicing occurs at this locus and four transcript variants encoding the same protein have been identified.




[View online »](#)


Synonyms: GS2; HsT18676; RAB27; RAM

Protein Families: Druggable Genome


Product images:



Gel: 8%SDS-PAGE
Lysate: 40 μ g
Lane: K562 cells
Primary antibody: TA350319 (RAB27A Antibody)
at dilution 1/400
Secondary antibody: Goat anti rabbit IgG at
1/8000 dilution
Exposure time: 30 seconds



Immunohistochemistry of paraffin-embedded
Human liver cancer tissue using TA350319
(RAB27A Antibody) at dilution 1/30 (Original
magnification: \times 200)



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using TA350319 (RAB27A Antibody) at dilution 1/30, treated with fusion protein. (Original magnification: $\times 200$)