

Product datasheet for **TA350471S**

Syntaxin 12 (STX12) Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1000-5000 WB positive control: Mouse brain tissue, hepg2, 293T and 231 cell
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Full length fusion protein
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	32 kDa
Gene Name:	syntaxin 12
Database Link:	NP_803173 Entrez Gene 23673 Human Q86Y82
Background:	Syntaxin 12, also known as STX12, STX13 or STX14, is a 276 amino acid single-pass membrane protein that contains one t-SNARE coiled-coil homology domain and belongs to the Syntaxin family. Syntaxin 12 regulates protein transport between late endosomes and the trans-Golgi network and interacts with ABC1 (ATP-binding cassette transporter A1), a protein that facilitates cellular release of choline-phospholipids and cholesterol to apoA-I (apolipoprotein A-I).
Synonyms:	STX13; STX14
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)

Protein Pathways: SNARE interactions in vesicular transport

Product images:



Gel: 8%SDS-PAGE
Lysate: 40 μ g
Lane 1-4: Mouse brain tissue
hepg2
293T and 231 cell
Primary antibody: [TA350471] (STX12 Antibody) at
dilution 1/800
Secondary antibody: Goat anti rabbit IgG at
1/8000 dilution
Exposure time: 5 seconds