

## Product datasheet for **TA350497S**


### THSD1 Rabbit Polyclonal Antibody

#### Product data:


|                       |   |
|-----------------------|---|
| Product Type:         | Primary Antibodies  |
| Applications:         | IHC   |
| Recommended Dilution: | IHC: 25-100<br>Positive control: Human thyroid cancer<br>Predicted cell location: Cytoplasm   |
| Reactivity:           | Human, Mouse  |
| Host:                 | Rabbit  |
| Isotype:              | IgG   |
| Clonality:            | Polyclonal  |
| Immunogen:            | Fusion protein of human THSD1   |
| Formulation:          | pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol  |
| Purification:         | Antigen affinity purification   |
| Conjugation:          | Unconjugated  |
| Storage:              | Store at -20°C as received.   |
| Stability:            | Stable for 12 months from date of receipt.  |
| Gene Name:            | thrombospondin type 1 domain containing 1   |
| Database Link:        | <a href="#">NP_954872</a><br><a href="#">Entrez Gene 56229 Mouse</a> <a href="#">Entrez Gene 55901 Human</a><br><a href="#">Q9NS62</a>  |
| Background:           | The protein encoded by this gene contains a type 1 thrombospondin domain, which is found in a number of proteins involved in the complement pathway, as well as in extracellular matrix proteins. Alternatively spliced transcript variants encoding different isoforms have been observed for this gene. |
| Synonyms:             | TMTSP; UNQ3010  |
| Protein Families:     | Transmembrane   |



[View online »](#)

**Product images:**

Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using [TA350497] (THSD1 Antibody) at dilution 1/25 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using [TA350497] (THSD1 Antibody) at dilution 1/25, treated with fusion protein. (Original magnification: ×200)