

Product datasheet for **TA350607S**

VPS53 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 200-1000 WB positive control: A172 cell IHC: 25-100 Positive control: Human lung cancer Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human VPS53
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	80 kDa
Gene Name:	VPS53, GARP complex subunit
Database Link:	NP_060759 Entrez Gene 68299 Mouse Entrez Gene 55275 Human Q5VIR6


Background: This gene encodes a protein with sequence similarity to the yeast Vps53p protein. Vps53p is involved in retrograde vesicle trafficking in late Golgi. May be involved in retrograde transport of early and late endosomes to the late Golgi. The GARP complex is required for the maintenance of the cycling of mannose 6-phosphate receptors between the TGN and endosomes, this cycling is necessary for proper lysosomal sorting of acid hydrolases such as CTSD.




[View online »](#)

Synonyms: HCCS1; hVps53L; PCH2E; pp13624


Product images:



Gel: 6%SDS-PAGE
Lysate: 40 µg
Lane: A172 cell
Primary antibody: [TA350607] (VPS53 Antibody)
at dilution 1/200
Secondary antibody: Goat anti rabbit IgG at
1/8000 dilution
Exposure time: 30 seconds



Immunohistochemistry of paraffin-embedded
Human lung cancer tissue using [TA350607]
(VPS53 Antibody) at dilution 1/20 (Original
magnification: ×200)



Immunohistochemistry of paraffin-embedded Human lung cancer tissue using [TA350607] (VPS53 Antibody) at dilution 1/20, treated with fusion protein. (Original magnification: ×200)