

Product datasheet for **TA350669S**

FOXC2 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 500-2000 WB positive control: Human breast infiltrative duct tissue
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide of human FOXC2
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	54 kDa
Gene Name:	forkhead box C2
Database Link:	NP_005242 Entrez Gene 14234 Mouse Entrez Gene 171356 Rat Entrez Gene 2303 Human Q99958
Background:	This gene belongs to the forkhead family of transcription factors which is characterized by a distinct DNA-binding forkhead domain. The specific function of this gene has not yet been determined; however, it may play a role in the development of mesenchymal tissues.
Synonyms:	FKHL14; LD; MFH-1; MFH1
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)

Product images:


Gel: 8%SDS-PAGE

Lysate: 40 µg

Lane: Human breast infiltrative duct tissue

Primary antibody: [TA350669] (FOXC2 Antibody)
at dilution 1/833

Secondary antibody: Goat anti rabbit IgG at
1/8000 dilution

Exposure time: 5 minutes