

Product datasheet for **TA350917S**

PIP5KI gamma (PIP5K1C) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 500-2000 WB positive control: A431 and hela cells IHC: 25-100 Positive control: Human liver cancer Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide of human PIP5K1C
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	73 kDa
Gene Name:	phosphatidylinositol-4-phosphate 5-kinase type 1 gamma
Database Link:	NP_036530 Entrez Gene 18717 Mouse Entrez Gene 314641 Rat Entrez Gene 23396 Human O60331

Background: This locus encodes a type I phosphatidylinositol 4-phosphate 5-kinase. The encoded protein catalyzes phosphorylation of phosphatidylinositol 4-phosphate, producing phosphatidylinositol 4,5-bisphosphate. This enzyme is found at synapses and has been found to play roles in endocytosis and cell migration. Mutations at this locus have been associated with lethal congenital contractural syndrome. Alternatively spliced transcript variants encoding different isoforms have been described.




[View online »](#)

Synonyms: LCCS3; PIP5K-GAMMA; PIP5K1-gamma; PIP5Kgamma


Protein Families: Druggable Genome

Protein Pathways: Endocytosis, Fc gamma R-mediated phagocytosis, Focal adhesion, Inositol phosphate metabolism, Metabolic pathways, Phosphatidylinositol signaling system, Regulation of actin cytoskeleton


Product images:




Gel: 10%SDS-PAGE
Lysate: 40 μ g
Lane 1-2: A431 cells
hela cells
Primary antibody: [TA350917] (PIP5K1C Antibody)
at dilution 1/500
Secondary antibody: Goat anti rabbit IgG at
1/8000 dilution
Exposure time: 2 minutes




Immunohistochemistry of paraffin-embedded
Human liver cancer tissue using [TA350917]
(PIP5K1C Antibody) at dilution 1/25 (Original
magnification: \times 200)



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using [TA350917] (PIP5K1C Antibody) at dilution 1/25, treated with synthetic peptide. (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human brain tissue using [TA350917] (PIP5K1C Antibody) at dilution 1/25 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human brain tissue using [TA350917] (PIP5K1C Antibody) at dilution 1/25, treated with synthetic peptide. (Original magnification: ×200)