

Product datasheet for **TA350947**

Artemin (ARTN) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 200-1000 WB positive control: 293T and SP20 cells, human fetal liver tissue
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide of human ARTN
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	23 kDa
Gene Name:	artemin
Database Link:	NP_476501 Entrez Gene 9048 Human Q5T4W7

Background: The protein encoded by this gene is a member of the glial cell line-derived neurotrophic factor (GDNF) family of ligands which are a group of ligands within the TGF-beta superfamily of signaling molecules. GDNFs are unique in having neurotrophic properties and have potential use for gene therapy in neurodegenerative disease. Artemin has been shown in culture to support the survival of a number of peripheral neuron populations and at least one population of dopaminergic CNS neurons. Its role in the PNS and CNS is further substantiated by its expression pattern in the proximity of these neurons.




[View online »](#)

Synonyms: artemin; ENOVIN; EVN; NBN; neublastin; neurotrophic factor; neurotrophic factor artemin; OTTHUMP00000009173; OTTHUMP00000009174

Protein Families: Druggable Genome, Secreted Protein

Product images:



Gel: 12%SDS-PAGE
Lysate: 50 µg
Lane 1-3: 293T cells
SP20 cells
human fetal liver tissue
Primary antibody: TA350947 (ARTN Antibody) at
dilution 1/350
Secondary antibody: Goat anti rabbit IgG at
1/8000 dilution
Exposure time: 2 minutes