

Product datasheet for TA351071

Froduct datasifeet for TASS 107

CENPE Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: IHC, WB

Recommended Dilution: WB: 200-1000

WB positive control: Human vermiform appendix tissue lysate

IHC: 50-100

Positive control: Human liver cancer Predicted cell location: Cytoplasm

Reactivity: Human, Mouse

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

Immunogen: Synthetic peptide of human CENPE

Formulation: pH7.4 PBS, 0.05% NaN3, 40% Glyceroln

Concentration: lot specific

Purification: Antigen affinity purification

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 316 kDa

Gene Name: centromere protein E

Database Link: NP 001804

Entrez Gene 1062 Human

Q02224

Background: Centrosome-associated protein E is a kinesin-like motor protein that accumulates in the G2

phase of the cell cycle. Unlike other centrosome-associated proteins, it is not present during

interphase and first appears at the centromere region of chromosomes during

prometaphase. CENPE is proposed to be one of the motors responsible for mammalian

chromosome movement and/or spindle elongation.



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn


Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com





Synonyms: CENP-E; KIF10; MCPH13; PPP1R61

Protein Families: Druggable Genome, Stem cell - Pluripotency

Product images:








Immunohistochemistry of paraffin-embedded Human liver cancer tissue using TA351071 (CENPE Antibody) at dilution 1/20 (Original magnification: ×200)





Immunohistochemistry of paraffin-embedded Human liver cancer tissue using TA351071 (CENPE Antibody) at dilution 1/20, treated with synthetic peptide. (Original magnification: ×200)