

Product datasheet for **TA351146S**

DRAK1 (STK17A) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 500-2000 WB positive control: Raji and 293T cells IHC: 25-100 Positive control: Human gastric cancer Predicted cell location: Cytoplasm
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide of human STK17A
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	47 kDa
Gene Name:	serine/threonine kinase 17a
Database Link:	NP_004751 Entrez Gene 9263 Human Q9UEE5
Background:	This gene is a member of the DAP kinase-related apoptosis-inducing protein kinase family and encodes an autophosphorylated nuclear protein with a protein kinase domain. The protein has apoptosis-inducing activity.
Synonyms:	DRAK1
Protein Families:	Druggable Genome, Protein Kinase




[View online »](#)


Product images:




Gel: 10%SDS-PAGE
Lysate: 40 μ g
Lane 1-2: Raji cells
293T cells
Primary antibody: [TA351146] (STK17A Antibody)
at dilution 1/900
Secondary antibody: Goat anti rabbit IgG at
1/8000 dilution
Exposure time: 1 minute




Immunohistochemistry of paraffin-embedded
Human gastric cancer tissue using [TA351146]
(STK17A Antibody) at dilution 1/25 (Original
magnification: \times 200)



Immunohistochemistry of paraffin-embedded Human gastric cancer tissue using [TA351146] (STK17A Antibody) at dilution 1/25, treated with synthetic peptide. (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human breast cancer tissue using [TA351146] (STK17A Antibody) at dilution 1/25 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human breast cancer tissue using [TA351146] (STK17A Antibody) at dilution 1/25, treated with synthetic peptide. (Original magnification: ×200)