

Product datasheet for **TA351229S**

GPR124 (ADGRA2) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 500-2000 WB positive control: Rat urinary bladder tissue lysate IHC: 50-200 Positive control: Human liver cancer Predicted cell location: Cell membrane
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide of human GPR124
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	143 kDa
Gene Name:	adhesion G protein-coupled receptor A2
Database Link:	NP_116166 Entrez Gene 78560 Mouse Entrez Gene 100363275 Rat Entrez Gene 25960 Human Q96PE1

Background: TEM5 (G-protein coupled receptor 124) and TEM7 (plexin domain containing 1). TEMs are associated with the cell surface membrane at low levels in normal human and mouse tissues. TEM5 is a seven-pass transmembrane receptor, whereas TEM1, TEM7 and TEM8 span the membrane once. TEM5 expression is elevated during tumor angiogenesis and neoangiogenesis. TEM7 is highly expressed in tumor endothelium and neurons. Therefore, TEM5 and TEM7 may be suitable targets for the development of antiangiogenic therapies.




[View online »](#)


Synonyms: GPR124; TEM5

Protein Families: Druggable Genome, Transmembrane


Product images:




Gel: 6%SDS-PAGE
Lysate: 40 µg
Lane: Rat urinary bladder tissue lysate
Primary antibody: [TA351229] (ADGRA2 Antibody) at dilution 1/400
Secondary antibody: Goat anti rabbit IgG at 1/5000 dilution
Exposure time: 5 minutes




Immunohistochemistry of paraffin-embedded Human liver cancer tissue using [TA351229] (ADGRA2 Antibody) at dilution 1/50 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using [TA351229] (ADGRA2 Antibody) at dilution 1/50, treated with synthetic peptide. (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human prostate cancer tissue using [TA351229] (ADGRA2 Antibody) at dilution 1/50 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human prostate cancer tissue using [TA351229] (ADGRA2 Antibody) at dilution 1/50, treated with synthetic peptide. (Original magnification: ×200)