

Product datasheet for **TA351260S**

PLAAT2 Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 200-1000 WB positive control: Human fetal liver tissue IHC: 100-300 Positive control: Human liver cancer Predicted cell location: Cytoplasm
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide of human HRASLS2
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	17 kDa
Gene Name:	HRAS like suppressor 2
Database Link:	NP_060348 Entrez Gene 54979 Human Q9NWW9
Background:	HRASLS2 (HRAS-like suppressor 2), also known as PLA1/2-2, is a 162 amino acid Cytoplasmic protein that belongs to the H-rev107 family. Expressed in liver, kidney, small intestine and colon, HRASLS2 exhibits PLA1/2 activity by catalyzing the calcium-independent hydrolysis of acyl groups in various phosphotidylcholines (PCs) and phosphatidylethanolamine (PE).
Synonyms:	2-2; PLA1
Protein Families:	Druggable Genome, Transmembrane




[View online »](#)

Product images:


Gel: 10%SDS-PAGE
Lysate: 40 μ g
Lane: Human fetal liver tissue
Primary antibody: [TA351260] (PLAAT2 Antibody) at dilution 1/362.5
Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution
Exposure time: 40 seconds




Immunohistochemistry of paraffin-embedded Human liver cancer tissue using [TA351260] (PLAAT2 Antibody) at dilution 1/40 (Original magnification: \times 200)



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using [TA351260] (PLAAT2 Antibody) at dilution 1/40, treated with synthetic peptide. (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using [TA351260] (PLAAT2 Antibody) at dilution 1/40 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using [TA351260] (PLAAT2 Antibody) at dilution 1/40, treated with synthetic peptide. (Original magnification: ×200)