

Product datasheet for TA351305

ITPR2 Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: IHC

Recommended Dilution: IHC: 50-200

Positive control: Human colon cancer Predicted cell location: Cytoplasm

Reactivity: Human
Host: Rabbit
Isotype: IgG

Clonality: Polyclonal

Immunogen: Synthetic peptide of human ITPR2

Formulation: pH7.4 PBS, 0.05% NaN3, 40% Glyceroln

Concentration: lot specific

Purification: Antigen affinity purification

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Gene Name: inositol 1,4,5-trisphosphate receptor type 2

Database Link: NP 002214

Entrez Gene 3709 Human

Q14571

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn





Background:

Inositol 1,4,5-triphosphate (IP3) functions as a second messenger for a myriad of extracellular stimuli including hormones, growth factors and neurotransmitters. Receptor tyrosine kinases indirectly increase the intracellular levels of IP3 through the activation of phospholipases such as phospholipase C (PLC), which convert phosphatidylinositol-4,5 bisphosphate into IP3 and diacylglycerol (DAG). The inositol 1,4,5-triphosphate receptor, IP3R, acts as an inositol triphosphate (IP3)-gated calcium release channel in a variety of cell types. Three IP3 receptor subtypes have been described and are designated IP3R-I, IP3R-II and IP3R-III. IP3R-I is the predominant IP3R subtype expressed in neuronal tissues and the central nervous system, but is also expressed at high levels in the liver.

Synonyms: CFAP48; IP3R2


Protein Families: Druggable Genome, Transmembrane

Protein Pathways: Alzheimer's disease, Calcium signaling pathway, Gap junction, GnRH signaling pathway, Long-

term depression, Long-term potentiation, Oocyte meiosis, Phosphatidylinositol signaling


system, Vascular smooth muscle contraction

Product images:



Immunohistochemistry of paraffin-embedded Human colon cancer tissue using TA351305 (ITPR2 Antibody) at dilution 1/40 (Original magnification: ×200)





Immunohistochemistry of paraffin-embedded Human colon cancer tissue using TA351305 (ITPR2 Antibody) at dilution 1/40, treated with synthetic peptide. (Original magnification: ×200)