

Product datasheet for TA351340

LARG (ARHGEF12) Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: IHC

Recommended Dilution: IHC: 50-200

Positive control: Human colon cancer Predicted cell location: Cytoplasm

Reactivity: Human, Mouse

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

Immunogen: Synthetic peptide of human ARHGEF12
Formulation: pH7.4 PBS, 0.05% NaN3, 40% Glyceroln

Concentration: lot specific

Purification: Antigen affinity purification

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Gene Name: Rho guanine nucleotide exchange factor 12

Database Link: NP 056128

Entrez Gene 69632 MouseEntrez Gene 23365 Human

Q9NZN5

Background: Rho GTPases play a fundamental role in numerous cellular processes that are initiated by

extracellular stimuli working through G protein-coupled receptors. The encoded protein may form a complex with G proteins and stimulate Rho-dependent signals. This protein has been

observed to form a myeloid/lymphoid fusion partner in acute myeloid leukemia. Two

transcript variants encoding different isoforms have been found for this gene.

Synonyms: LARG; PRO2792

Protein Pathways: Axon guidance, Regulation of actin cytoskeleton, Vascular smooth muscle contraction




OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn


Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



Product images:




Immunohistochemistry of paraffin-embedded Human colon cancer tissue using TA351340 (ARHGEF12 Antibody) at dilution 1/30 (Original magnification: ×200)




Immunohistochemistry of paraffin-embedded Human colon cancer tissue using TA351340 (ARHGEF12 Antibody) at dilution 1/30, treated with synthetic peptide. (Original magnification: ×200)





Immunohistochemistry of paraffin-embedded Human breast cancer tissue using TA351340 (ARHGEF12 Antibody) at dilution 1/30 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human breast cancer tissue using TA351340 (ARHGEF12 Antibody) at dilution 1/30, treated with synthetic peptide. (Original magnification: ×200)