

Product datasheet for **TA351741S**

CCS Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 500-2000 WB positive control: Human fetal liver tissue IHC: 25-100 Positive control: Human liver cancer Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide of human CCS
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	29 kDa
Gene Name:	copper chaperone for superoxide dismutase
Database Link:	NP_005116 Entrez Gene 12460 MouseEntrez Gene 84485 RatEntrez Gene 9973 Human O14618
Background:	Copper chaperone for superoxide dismutase specifically delivers Cu to copper/zinc superoxide dismutase and may activate copper/zinc superoxide dismutase through direct insertion of the Cu cofactor.
Synonyms:	MGC138260
Protein Families:	Druggable Genome




[View online »](#)

Protein Pathways: Amyotrophic lateral sclerosis (ALS)


Product images:



Gel: 10%SDS-PAGE
Lysate: 40 µg
Lane: Human fetal liver tissue
Primary antibody: [TA351741] (CCS Antibody) at dilution 1/300
Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution
Exposure time: 40 seconds



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using [TA351741] (CCS Antibody) at dilution 1/35 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using [TA351741] (CCS Antibody) at dilution 1/35, treated with synthetic peptide. (Original magnification: $\times 200$)