

Product datasheet for **TA351878S**

Thromboxane A2 receptor (TBXA2R) Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 500-2000 WB positive control: Raji and Jurkat cells IHC: 30-150 Positive control: Human cervical cancer Predicted cell location: Cell membrane
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide of human TBXA2R
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	37 kDa
Gene Name:	thromboxane A2 receptor
Database Link:	NP_963998 Entrez Gene 6915 Human P21731
Background:	This gene encodes a member of the G protein-coupled receptor family. The protein interacts with thromboxane A2 to induce platelet aggregation and regulate hemostasis. A mutation in this gene results in a bleeding disorder. Multiple transcript variants encoding different isoforms have been found for this gene.
Synonyms:	BDPLT13; TXA2-R
Protein Families:	Druggable Genome, GPCR, Transmembrane




[View online »](#)

Protein Pathways: Calcium signaling pathway, Neuroactive ligand-receptor interaction


Product images:




Gel: 8%SDS-PAGE
Lysate: 40 μ g
Lane 1-2: Raji cells
Jurkat cells
Primary antibody: [TA351878] (TBXA2R Antibody)
at dilution 1/200
Secondary antibody: Goat anti rabbit IgG at
1/8000 dilution
Exposure time: 40 seconds




Immunohistochemistry of paraffin-embedded
Human cervical cancer tissue using [TA351878]
(TBXA2R Antibody) at dilution 1/40 (Original
magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human cervical cancer tissue using [TA351878] (TBXA2R Antibody) at dilution 1/40, treated with synthetic peptide. (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using [TA351878] (TBXA2R Antibody) at dilution 1/40 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using [TA351878] (TBXA2R Antibody) at dilution 1/40, treated with synthetic peptide. (Original magnification: ×200)